RSC Advances



This is an *Accepted Manuscript*, which has been through the Royal Society of Chemistry peer review process and has been accepted for publication.

Accepted Manuscripts are published online shortly after acceptance, before technical editing, formatting and proof reading. Using this free service, authors can make their results available to the community, in citable form, before we publish the edited article. This Accepted Manuscript will be replaced by the edited, formatted and paginated article as soon as this is available.

You can find more information about *Accepted Manuscripts* in the **Information for Authors**.

Please note that technical editing may introduce minor changes to the text and/or graphics, which may alter content. The journal's standard <u>Terms & Conditions</u> and the <u>Ethical guidelines</u> still apply. In no event shall the Royal Society of Chemistry be held responsible for any errors or omissions in this *Accepted Manuscript* or any consequences arising from the use of any information it contains.



www.rsc.org/advances

Cite this: DOI: 10.1039/c0xx00000x

www.rsc.org/xxxxx

Hydrazine Gas

ARTICLE TYPE Synthesis of Dihydroquinoline Based Merocyanine as 'Naked Eye' and 'Fluorogenic' sensor for Hydrazine Hydrate in Aqueous Medium and

K.Vijay,^a C. Nandi,^b and Shriniwas D. Samant^{*a}

s Received (in XXX, XXX) Xth XXXXXXXX 20XX, Accepted Xth XXXXXXXX 20XX DOI: 10.1039/b000000x

1,2,2,4-Tetramethyldihydroquinoline (TMDQ) is formylated at the 6-position and the resulting quinoline aldehyde is condensed with 3-methyl-1-phenyl-1H-pyrazol-5(4H)-one to obtain a merocyanine 6pyrazolinyl TMDQ. The merocyanine is a selective and sensitive sensor for hydrazine hydrate in aqueous

10 ethanolic medium and also for hydrazine gas. The merocyanine detects hydrazine hydrate selectively in less than 5s over other amines, anions, and metals. The detection can be done simply by naked eye and quantification can be done by absorbance / fluorescence spectroscopy. Hydrazine hydrate as well as hydrazine gas can be detected by TLC plate technique.

Introduction

- 15 Hydrazine is an important industrial chemical used widely in rockets and missiles as fuel¹, as a reducing agent², and as an antioxidant in nuclear and electrical power plants¹. It has vital function in the pharmaceutical and chemical industries as catalyst, dyes in textiles, and inhibitors³. Despite of many
- 20 applications, hydrazine hydrate is extremely toxic. Hydrazine is most likely a carcinogen and its threshold limit value (TLV) is 10 ppb⁴. It can induce damage of the central nervous system (CNS), nose irritation, and temporary blindness⁵. Hence, designing a reliable and rapid detector that can provide on-site and real-time 25 recognition of hydrazine is very desirable.
- There are several reports on sensing of hydrazine⁶⁻¹³. But very few reports of 'Naked Eye' and 'Fluorescent' probes for hydrazine hydrate can be found in the literature. Zhao reported a hydrazine sensor which contains expensive gold metal¹¹. Further,
- 30 the sensor does not show any fluorescence properties. Very recently, a hydrazine sensor which shows fluorescence was reported.¹³ Similar molecule was reported by Cui almost at same time.¹⁴ However a probe which can be used either 'Naked Eye', or in 'UV-vis absorption', as well as in 'Fluorescence technique' 35 is definitely a superior one.
- Cyanines were initially used as developers in photography. Because of their conjugated structure and useful properties, people started using cyanines for diverse applications¹⁵. Merocyanines are a class of cyanines which have various
- ⁴⁰ applications.¹⁶ Merocyanines can be used as probes and markers in chemical analysis, biology and medical applications, nonlinear optics, and organic semiconductors¹⁶. Many cyanines contain quinoline moiety¹⁷
- ^a Department of Chemistry, Institute of Chemical Technology, Matunga, 45 Mumbai-400019, India. E-mail: samantsd@yahoo.com.
 - ^b R & D, NOCIL limited, C 37, TTC Industrial Area, Navi Mumbai-400705, India.

There are various detection methods for sensing toxic chemicals¹⁸, such as naked eye detection, spectrophotometry¹⁹, ⁵⁰ spectrofluorimetry²⁰, electrochemistry²¹, GC²², and HPLC²³. Each method has its own pros and cons with respect to portability, analysis time, sensitivity, and economy. However, a sensor which can be used in multiple detection methods can overcome these drawbacks and thus is advantageous. For example, 'naked

55 eye' sensors particularly show potential because the alteration in colour can easily be observed without any equipment and hence are more economic. On the other hand, UV instrument can permit the quantification by the measurement of the absorption intensity. Further, if a sensor shows fluorescence, it is an added benefit. For

- 60 example, the analyte may be determined with more accuracy in a spectrofluorimeter than in an UV instrument. In addition, fluorescence can be used in biological applications like cell imaging.
- Merocyanines have a potential as sensors and many of them 65 contain a quinoline nucleus. This prompted us to synthesize a merocyanine dye with dihydroquinoline joined to a heterocyclic ring through a conjugation extending linker. In this paper we would like to report the synthesis of (E)-3-methyl-1-phenyl-4-((1,2,2,4-tetramethyl-1,2-dihydroquinolin-6-yl)methylene)-1H-
- 70 pyrazol-5(4H)-one which belongs to the class of merocyanines and is a conjugate of 1,2,2,4-tetramethyldihydroquinoline and pyrazolone. The merocyanine can effectively and competently sense hydrazine hydrate. It can act as a rapid and selective nakedeye, fluorescent, as well as absorbance sensor for hydrazine 75 hydrate in an aqueous alcoholic medium. It is easy to synthesize, eco-friendly, sensitive, and selective. The sensing does not involve any complicated buffer-making procedure, and can be done simply in an ethanol: water system. To the best of our knowledge, this is the first report of a merocyanine acting as a 80 hydrazine sensor in an aqueous alcoholic medium and also for hydrazine gas.

This journal is © The Royal Society of Chemistry [2014]

Results and discussion

- 1,2,2,4-tetramethyl-1,2-dihydroquinoline(1) was prepared by the condensation of aniline with acetone²⁴. 1 was N-methylated by methyl iodide in DMF in the presence of potassium carbonate to ⁵ obtain 2. 2 was formylated using DMF/POCl₃ according to the reported procedure to 1,2,2,4-tetramethyl-1,2-dihydroquinoline-6-carbaldehyde (3)²⁴ Phenyl hydrazine was condensed with ethyl acetoacetate to obtain 3-methyl-1-phenyl-1H-pyrazol-5(4H)-ones (6) as per the reported procedure²⁵. 3 and 6 were condensed in the ¹⁰ presence of TEA in ethanol to obtain **7** in 82% yield (scheme-1).
- The structure of 7 was confirmed by ¹H NMR and ¹³C NMR as well as by ESI-MS. It was an intensely orange coloured powder



Scheme 1 Synthesis of (E)-3-methyl-1-phenyl-4-((1,2,2,4-tetramethyl-1,2-dihydroquinolin-6-yl)methylene)-1H-pyrazol-5(4H)-one.

Water is the best solvent for any sensor. 7 was insoluble in water, but soluble completely in ethanol-water (8:2 v/v). Hence, all ³⁰ studies were done in ethanol-water (8:2 v/v). The UV-vis spectrum showed λ_{max} at 503nm. When increasing concentration of hydrazine hydrate was put in to a fixed concentration of 7, the absorbance of the dye decreased. Hydrazine hydrate at different concentrations was added to 16.16 μ M (concentration which

- ³⁵ obeyed Beer-Lambert's law) of probe-7, (Fig. 1).This showed that with increasing concentration of hydrazine hydrate the red color of dye decreased linearly and finally vanished. This visual response was instantaneous (<5 seconds). Thus, 7 is a 'nakedeye' colorimetric sensor. Moreover, the probe was not fluorescent initially, but after addition of hydrazine hydrate, it showed
- ⁴⁰ initially, but after addition of hydrazine hydrate, it showed fluorescence at 432 nm.



Fig. 1 Visual colour change in 7 (16.16 μ M, in ethanol: water (8:2)) when hydrazine hydrate is added at different concentrations (0, 2, 3, and 4 mM, a-d respectively).

The colour of the merocyanine is due to extended conjugation from quinoline ring to pyrazoline ring through the exo-cyclic double bond on the pyrazoline ring. Hydrazine hydrate may attack the linker double bond that is conjugated to the carbonyl ⁵⁵ group in pyrazolone and cleave the molecule to form hydrazone 8 and pyrazolone 6 (scheme 2). This leads to a change in the absorbance of the probe. To confirm this, the probe-7 was mixed with 5 eq of hydrazine hydrate in ethanol and stirred for 1 h at room temperature and the reaction was monitored by TLC which ⁶⁰ showed the formation of a new compound. This compound was matching with hydrazone of quinoline aldehyde 3. The reaction was worked up and purified by flash column chromatography_to get an off white to light yellow coloured product. The structure was further confirmed by ¹H, ¹³C NMR.



Scheme 2 Proposed reaction of hydrazine with probe-7

70 UV-Vis study

The probe sensed hydrazine hydrate very rapidly within few seconds and a visible response was seen. In order to follow the reaction between probe-7 and hydrazine hydrate, the absorbance of 7 at 503 nm was recorded with respect to time after mixing ⁷⁵ hydrazine hydrate. Using 160 x 10^{-5} M hydrazine hydrate (concentration where the absorbance is almost negligible) and 16.16 μ M of 7, absorbance was measured with respect to time (Fig.2). The initial response was very rapid. The reaction was complete after 90 min.



Fig. 2 Plot of absorbance at 503 nm with respect to time of a system containing 160×10^{-5} M hydrazine hydrate and 16.16μ M of 7.

⁹⁵ UV–Vis spectrum of **7** was recorded with increasing concentration of hydrazine hydrate (Fig.3a and 3b). Addition of hydrazine to the solution of **7** led to a decrease in the absorbance, which vanished after addition of about 99 equiv (160×10^{-5} M) of hydrazine hydrate. All the spectra were recorded 90 min after the ¹⁰⁰ addition of hydrazine hydrate. In the range of 2 x 10^{-5} to 180 x 10^{-5} M of hydrazine hydrate for 16.16 µM of **7**, the change in absorption was inversely proportional to the concentration of hydrazine. Thus, **7** acted as an efficient 'turn-off absorbance' probe for the selective detection of hydrazine hydrate in an ¹⁰⁵ aqueous alcoholic medium. The experiments were repeated for three times to check the reproducibility of the results.

A graph of extinction ratio $(A_{503}^0 - A_{503}^0)/A_{503}^0$ versus hydrazine hydrate concentration was plotted (Fig.3c and 3d).This showed good sensitivity of the probe upon addition of hydrazine hydrate.

45

(ESI[†]).

The extinction ratio showed a good linear relationship with concentration of hydrazine hydrate between 1.23 to 61.9 equivalents (i.e., 2 to 100 x 10^{-5} M) (Co-efficient of determination (R²) given in the graph). The detection limit²⁶ of ⁵ probe-7 for hydrazine hydrate was estimated to be of 16×10^{-7} M







Fig.4(a) UV response of probe-7 (16.16 μM) to 50 equivalents of various metals and hydrazine. (b) UV response of probe-7 (16.16 μM) to 50 equivalents of various anions and hydrazine. (C) UV response of probe-7 (16.16 μM) to 50 equivalents of various neutral molecules.

50

Selectivity of 7

Some metals may form complex with the probe and interfere with the sensing process. In fact anions and amines also may interfere. Hence the selectivity of the probe towards hydrazine

⁵⁵ hydrate was studied with respect to various common metals, anions, and neutral molecules. For 1 equivalent of probe-7, 50 equivalents of interfering anions/cations/neutral molecules were added separately and the absorption spectra was recorded (ESI[†]). All these spectra were recorded 90 minutes after the ⁶⁰ addition of anions/cations/neutral molecules. The graph showing the extinction ratio at 503 nm was plotted (Fig.4a, 4b, & 4c). Hydrazine hydrate (50 equivalents / 80 x 10⁻⁵ M) only is found to affect the absorption of 7 to a significant extent. Thus probe-7 showed excellent selectivity towards hydrazine ⁶⁵ hydrate.

Fluorescence study

An interesting feature of probe-7 was that after addition of hydrazine hydrate is UV absorption was '*turned off*'. When the resulting solution was kept in an UV chamber at 254 nm, ⁵ hydrazine hydrate '*turned on*' the emission feature of the probe. Hence the fluorescence study of the probe was carried out. The fluorescence spectra of probe-7 and 7 in the presence of increasing concentrations of hydrazine hydrate are given in Fig.5a. The action of hydrazine with probe-7 turns-on the ¹⁰ emission feature and shows peak at 432 nm. Fluorescence

- spectral changes were estimated after 90 minutes quantitatively with increasing concentrations of hydrazine (20–400 x 10^{-5} M) in ethanol: water system (8:2) (Fig.5b). The fluorescence intensity reached a point of negligible deviation upon the
- ¹⁵ addition of 99 equiv. (160.0 x 10^{-5} M) of hydrazine. The ratiometric changes under fluorescence were same as those observed in UV-vis studies. A linear relationship was found between the emission intensity at 432 nm and the hydrazine concentration over 20–140 x 10^{-5} M range. The fluorescence
- ²⁰ response of the probe with respect to other common interfering anions, cations, and related amino molecules also was checked. This study showed very high selectivity of the probe towards hydrazine hydrate among others. (ESI⁺)



Fig.5(a) Fluorescence spectra of probe-7 (16.16 μ M) with increasing concentrations of hydrazine hydrate. (b) Variation of Fl intensity of probe-7 (16.16 μ M) at 432 nm with increasing concentrations of hydrazine hydrate; error bars for extinction ratio are shown using standard deviations from the mean.

NMR study:

⁶⁰ ¹HNMR and ¹³CNMR were used for characterizing the product. The ¹HNMR of the probe-7 and that of product-8 on reaction of 7 with hydrazine hydrate were recorded. From the spectra it appears that the pyrazolone part was entirely removed in **8**. (Fig.6). ¹HNMR of product-8 shows loss of methyl ($\delta = 2.33$) and phenyl groups of pyrazolone unit. The proton of probe-7 at 8.47 (=CH-C) is shifted to 8.58 (N=CH-C) in the product-8. Further when the ¹HNMR spectrum of **8** was compared with that of 3, the only difference was found on the peak position of aldehydic H of 3 ($\delta = 9.8$) and the proton of imine ($\delta = 8.58$).



Fig. 6 1 H NMR spectra in CDCl₃ of probe-7 and Product-8

90

105

TLC plate detection techniques:

Detection of hydrazine hydrate by TLC plates:

Probe-7 was dissolved in ethyl acetate to get a 5 mM solution. Thin layer chromatography (TLC) plates were dipped ⁹⁵ into it and kept for 30 min at room temperature for drying. These plates were immersed into solutions of different concentrations of hydrazine hydrate in ethanol: water (8:2) medium and then exposed to air to get rid of the solvent. A minimum of 1 mM of hydrazine hydrate was easily ¹⁰⁰ recognizable by naked eye in this TLC plate technique. This is similar as observed in naked eye detection in aqueous ethanolic medium in fig. 1. Hence the detection limit of hydrazine hydrate by TLC plate technique is 1 mM. In Fig. 7, the complete colour change from orange to colourless is shown.



Detection of hydrazine gas:

The probe could rapidly and effectively sense hydrazine gas. ¹⁰ This experiment was carried out using a TLC plate coated with probe-7. A 5 mM solution of the probe was prepared in ethyl acetate. A TLC plate was dipped into it and dried at room temperature for 30 min. The TLC plate was kept in a vial. Excess of hydrazine gas was passed into the vial for 5 sec. The ¹⁵ colour of the TLC plate changed from orange to colourless within 5 sec. (Fig.8).



Fig.8 TLC plate before hydrazine gas passing, and after excess hydrazine gas passing for 5 seconds.

DFT study:

20

25

In order to throw more light on the actual chromo/fluorophore, DFT studies were done on 7 and the possible reduced product-8. DFT study was done using Gaussian 03 program²⁷. We ³⁰ performed B3LYP exchange correlation functional under 6-

- 311G+, p, d basis set and calculated HOMO and LUMO levels. The optimized geometries and calculated electron distributions in the FMOs of probe-7 and product-8 are shown in Fig. 9. The HOMO-LUMO orbital distribution is spread from the quinoline ³⁵ moiety to the pyrazolone moiety. There is a planarity over the
- ³⁵ molecty to the pyrazolone molecy. There is a planality over the major chromophore extending from quinoline to pyrazolone unit. The HOMO-LUMO energy gap in the probe-7 is less than that of the product-8 because of extended conjugation. This observation is matching with the fact that the probe is coloured ⁴⁰ and the product is not.

Conclusion

50

A merocyanine probe is described that provides very good sensitivity for detecting hydrazine without needing any advanced, complex readout equipment or complex buffer ⁴⁵ making procedure. The synthesis of the probe is very easy. The detection limit of the probe for hydrazine hydrate is 1.6 μ M. The gaseous hydrazine is detected by the probe within 5 seconds. The TLC plate technique provides a more practical and sensitive detection method.



70 Fig. 9 HOMO-LUMO levels and optimized structures of probe-7 and its reduced product from Gaussian 03.

Acknowledgement

The authors are thankful to IIRBS, Kerala, India and GJUST, Punjab, India for providing NMR spectra. KV is thankful for 75 CSIR, New Delhi, India for a fellowship.

Notes and references

† Electronic Supplementary Information (ESI) available. See DOI: 10.1039/b000000x/

80 Reference List

(1) Zelnick, S. D.; Mattie, D. R.; Stepanaik, P. C. Aviation, Space and environmental Medicine 2003, 74, 1285-1291.

- (2) Ragnarsson, U. Chem.Soc.Rev. 2001, 30, 205-213.
- 85 (3) Narayanan, S. S.; Scolz, F. *Electroanalysis* 1999, 11, 465-469.

(4) Documentation of the Threshold Limit Values and Biological Exposure Indices. American Conference of Governmental Industrial Hygenists. 1999.

90 Ref Type: Conference Proceeding (5) Wang G: Zhang C: He X: I

(5) Wang, G.; Zhang, C.; He, X.; Li, Z.; Zhang, X.; Wang, L.; Bang, B. *Electrochim.Acta* 2010, *55*, 7204-7210.
(6) Goswami, S.; paul, S.; Manna, A. *RSC Adv.* 2013, *3*, 18872-

18877.
(7) Lee, M. H.; Yoon, B.; Kim, J. S.; Sessler, J. L. *Chem.Sci.* 2013, *4*, 4121-4126.

(8) Li, K.; Xu, H.-R.; Yu, K.-K.; Hou, J.-T.; Yu, X.-Q. Anal.methods **2013**, *5*, 2653-2656.

(9) Virji, S.; Kaner, R. B.; Weiller, B. H. *Chem.Mater.* 2005, *17*, 1256-1260.

(10) Wang, J.; Chen, L. Anal. Chem. 1995, 67, 3824-3827.

- (11) Zhao, Z.; Zhang, G.; Gao, Y.; Yang, X.; Li, Y. *Chem. Commun.* **2011**, *47*, 12816-12818.
- (12) Lin, Y. D.; Chow, T. J. RSC Adv. 2013, 3, 17924-17929.

(13) Ramakrishnam Raju, M. V.; Chandra Prakash, E.; Chang, H. C.; Lin, H. C. *Dyes and Pigments* **2014**, *103*, 9-20.

- (14) Cui, L.; Peng, Z.; Ji, C.; Huang, J.; Huang, D.; Ma, J.; Zhang, S.; Qian, X.; Xu, Y. Chem. Commun. **2014**, *50*, 1485-1487.
- 5 (15) Mishra, A.; Bahera, R. K.; Bahera, P. K.; Mishra, B. K.; Bahera, G. B. *Chem.Rev.* 2000, 100, 1973-2011.
 (16) Kulinich, A. V.; Ishchenko, A. A. *Russian Chemical Reviews* 2009, 78, 141-164.
 (17) Deligeorgiev, T. G.; Kaloyanova, S.; Vaquero, J. J. *Recent*
- Deligeorgiev, I. G.; Kaloyanova, S.; Vaquero, J. J. Recent
 Patents on Materials Science 2009, 2, 1-26.
 - (18) Abbaspour, A.; Mirahmadi, E.; Khajehzadeh, A. *Anal.methods* **2010**, *2*, 349-353.
 - (19) Afkhami, A.; Bahram, M. *Talanta* **2006**, *68*, 1148-1155.
 - (20) Ensafi, A. A.; Rezaei, B. *Talanta* **1998**, *47*, 645-649.
 - (21) Batchelor-McAuley, C.; Banks, C. E.; Simm, A. O.; Jones, T. G.; Compton, R. G. *Analyst* **2006**, *131*, 106-110.
 (22) Sun, M.; Bai, L.; Liu, D. *J.Pharm.Biomed.Anal.* **2006**, *49*, 529-533.
- (23) Bhutani, H.; Singh, S.; Vir, S.; Bhutani, K. K.; Kumar, R.;
 Chakraborti, A. K.; Jindal, K. C. *J.Pharm.Biomed.Anal.* 2007, *43*,
 - (24) Tian, K.; Hu, D.; Hu, R.; Wang, S.; Li, S.; Li, Y.; Yang, G.
 - (24) Tian, K.; Hu, D.; Hu, K.; Wang, S.; Li, S.; Li, Y.; Yang, G. *Chem.Commun.* **2011**, *47*, 10052-10054.
 - (25) Nayak, M.; Batchu, H.; Batra, S. Tetrahedron Lett. 2012, 53, 4206-4208.
 - (26) Long, G. L.; Winefordner, J. D. Anal.Chem. 1983, 55, 712A-724A.

(27) Frisch, M. j. Gaussian 03[Revision B.03]. 2003. Pittsburgh, PA. Ref Type: Computer Program

30

25

15