

# Journal of Materials Chemistry B

Materials for biology and medicine

rsc.li/materials-b



Themed issue: Emerging Investigators 2020

ISSN 2050-750X

## PAPER

Ann-Christin Pöppler *et al.*  
 $^{14}\text{N}$ - $^1\text{H}$  HMQC solid-state NMR as a powerful tool to study  
amorphous formulations – an exemplary study of paclitaxel  
loaded polymer micelles

## PAPER

[View Article Online](#)  
[View Journal](#) | [View Issue](#)Cite this: *J. Mater. Chem. B*, 2020, 8, 6827 **$^{14}\text{N}$ – $^1\text{H}$  HMQC solid-state NMR as a powerful tool to study amorphous formulations – an exemplary study of paclitaxel loaded polymer micelles†**Marvin Grüne,<sup>a</sup> Robert Luxenhofer,<sup>b</sup> Dinu Iuga,<sup>c</sup> Steven P. Brown<sup>c</sup> and Ann-Christin Pöppler<sup>id</sup> \*<sup>a</sup>

Amorphous drug–polymer formulations are complex materials and often challenging to characterize, even more so if the small molecule component itself is increasingly complex. In this work, we present  $^{14}\text{N}$ – $^1\text{H}$  HMQC magic-angle spinning (MAS) NMR experiments in the solid state as a promising tool to study amorphous formulations. Poly(2-oxazoline) based polymer micelles loaded with different amounts of the cancer drug paclitaxel serve to highlight the possibilities offered by these experiments: while the dense core of these polymeric micelles prevents NMR spectroscopic analysis in solution and the very similar  $^{15}\text{N}$  chemical shifts hamper a solid-state NMR characterization based on this nucleus,  $^{14}\text{N}$  is a very versatile alternative.  $^{14}\text{N}$ – $^1\text{H}$  HMQC experiments yield well-separated signals, which are spread over a large ppm range, and provide information on the symmetry of the nitrogen environment and probe  $^{14}\text{N}$ – $^1\text{H}$  through-space proximities. In this way, the overall complexity can be narrowed down to specific N-containing environments. The results from the experiments presented here represent a valuable puzzle piece, which helps to improve the structural understanding of drug–polymer formulations. It can be straightforwardly combined with complementary NMR spectroscopic experiments and other analytical techniques.


Received 5th March 2020,  
Accepted 28th May 2020

DOI: 10.1039/d0tb00614a

[rsc.li/materials-b](http://rsc.li/materials-b)**Introduction**

Paclitaxel (PTX, Scheme 1) is an effective anti-cancer drug for a wide range of tumours, but it exhibits a very low aqueous solubility of  $0.4\ \mu\text{g mL}^{-1}$ .<sup>1,2</sup> Therefore, a variety of different formulations for PTX has been developed.<sup>3</sup> This includes the protein-based nanoparticle Abraxane<sup>®</sup>, the polymer-conjugate Opaxio<sup>™</sup>, Genexol-PM<sup>®</sup> and NK105.<sup>4,5</sup> The latter two contain PTX in polymer micelles formed of polyethylene glycol–poly(lactic acid) (Genexol-PM<sup>®</sup>) and polyethylene glycol–poly(amino acid) (NK105) copolymers. In general, the development of formulations and their clinical trials require very long times with highly uncertain outcomes. The above mentioned formulations have undergone multiple phase II/III clinical trials over the past decades, *e.g.* a recently published phase III clinical trial of NK105,<sup>6</sup> and only Abraxane<sup>®</sup> has been successfully marketed with FDA approval since 2007. Furthermore, they are all

characterized by a relatively low PTX loading not exceeding 25 wt%. In this context, Luxenhofer *et al.* reported promising preclinical data for polymer micelles comprising a poly(2-oxazoline) (POx) based triblock copolymer loaded with up to 50 wt% of PTX.<sup>7,8</sup> Compared to the clinically approved PTX



**Scheme 1** Structural formula of the components used in this study: the amphiphilic block copolymer poly(2-methyl-2-oxazoline)-block-poly(2-*n*-butyl-2-oxazoline)-block-poly(2-methyl-2-oxazoline) (POL) encapsulates paclitaxel (PTX) by self-assembly into polymeric micelles (schematic drawing on the right).

<sup>a</sup> Institute of Organic Chemistry, University of Würzburg, Am Hubland, 97074 Würzburg, Germany. E-mail: ann-christin.poeppler@uni-wuerzburg.de<sup>b</sup> Lehrstuhl für Chemische Technologie der Materialsynthese, University of Würzburg, Röntgenring 11, 97070 Würzburg, Germany<sup>c</sup> Department of Physics, University of Warwick, Coventry, CV4 7AL, UK

† Electronic supplementary information (ESI) available. See DOI: 10.1039/d0tb00614a



formulations, their formulations showed higher maximum tolerated doses, elevated drug exposure to tumour tissues and prolonged survival of mice bearing A2780 human ovarian tumours. Schulz *et al.*,<sup>9</sup> Jaksch *et al.*<sup>10</sup> and Sochor *et al.*<sup>11</sup> studied their micellar morphology for different drug loadings (PTX and curcumin) using dynamic light scattering, atomic force microscopy, (cryogenic) transmission electron microscopy, and small-angle neutron scattering (SANS). They found that the pure polymer self-assembles into wormlike and spherical micelles in aqueous solution, while the incorporation of PTX led to the exclusive formation of spherical particles. Interestingly, SANS data revealed small, PTX rich domains partly submerged within the micellar core. Apart from such morphological aspects and capability to design improved drug delivery systems, it is very important to also understand the complex structural arrangements in such amorphous drug-polymer formulations on the molecular level. NMR probes the local environment of molecules and is thus an excellent tool to report on various possible interactions. This doesn't necessarily require crystalline or well-ordered structures, and also gives excellent results for disordered as well as amorphous samples.<sup>12–14</sup>

In particular, solid-state NMR has been established as a powerful tool for structure elucidation in pharmaceutical contexts such as identification of polymorphs or investigation of amorphous solid dispersions.<sup>15–17</sup> The sometimes strongly reduced mobility of encapsulated drugs within drug formulations can hinder the characterization by NMR in solution due to strong broadening of the signals, which makes solid-state NMR the corresponding method of choice. Callari *et al.* as well as Pöppler *et al.* recently showcased how solid-state NMR at moderate to fast Magic Angle Spinning (MAS) helps to obtain loading dependent structural insights into micellar formulations with an assumed core-shell structure.<sup>18,19</sup> Both groups found that increasing the loading of two different types of polymeric micelles did not just affect the micellar core but interestingly also affected the surrounding shell, which served as a basis to explain their physicochemical and biological properties such as reduced cellular uptake and inferior dissolution rates. Hirsh *et al.* have recently also shown how <sup>1</sup>H MAS measurements enable the rapid characterization of general pharmaceutical dosage forms.<sup>20</sup> The enormous progress made with respect to hardware development, *e.g.* MAS frequencies of up to 130 kHz being available now,<sup>21</sup> also enables high-resolution detection of protons in the solid state with longer coherence lifetimes.<sup>22–24</sup> The latter are essential prerequisites for sophisticated homonuclear as well as heteronuclear 2D experiments, which are exciting tools to learn more about the intermolecular proximities in a large range of materials.

However, for the solid-state NMR spectroscopic characterization of PTX and its formulations as discussed above, a series of complexities arises: PTX alone with its 51 protons and thus a variety of strongly overlapping chemical environments already poses a challenge. In addition, crystalline PTX contains two individual molecules in the asymmetric unit ( $Z' = 2$ ),<sup>25</sup> which leads to the doubling of the expected resonances. Consequentially, to date, only <sup>13</sup>C solid-state NMR data with assignment based on

comparison with structurally related fragments and calculations has been published,<sup>26,27</sup> while <sup>1</sup>H solid-state NMR data of PTX is still missing. Taking the polymer component and potential drug-polymer interactions into account further complicates the spectral evaluation and thus requires narrowing down the search space, *e.g.* by focusing on specific functionalities. For example, spectral simplification can be achieved by making use of NMR active heteronuclei.

The frequent involvement of nitrogen atoms in hydrogen bonding in pharmaceuticals makes this nucleus a promising starting point for a detailed investigation of PTX containing formulations through 2D nitrogen-proton correlation experiments, especially since PTX has only one nitrogen environment per molecule. There are two NMR-active nitrogen isotopes with low magnetogyric ratios: <sup>15</sup>N ( $I = 1/2$ ) possesses a low natural abundance of only 0.37%, whereas <sup>14</sup>N has a natural abundance of 99.6% but exhibits quadrupole broadening due to its nuclear spin quantum number  $I = 1$ . The substantial disadvantage of the first-order quadrupole interaction, causing a broadening in the range of several MHz, can be avoided by setting the indirect <sup>14</sup>N spectral width of the 2D correlation equal to the spinning frequency. Tatton *et al.* successfully applied this in the form of a <sup>14</sup>N-<sup>1</sup>H heteronuclear multiple-quantum correlation (HMQC) experiment to study crystalline model compounds,<sup>28</sup> hydrogen bonding in co-crystals and attempting the first transfer to amorphous solid dispersions.<sup>29</sup> So far, the majority of <sup>14</sup>N-<sup>1</sup>H HMQC experiments have been performed for highly ordered, crystalline compounds.<sup>28,30–39</sup> In contrast, disordered systems with significant differential dynamics are more challenging systems and there are thus very seldom reports on employing the <sup>14</sup>N-<sup>1</sup>H HMQC experiment for their characterization. However, the characterization of disordered and amorphous materials could substantially benefit from the additional <sup>14</sup>N second-order isotropic quadrupolar shifts (see the ESI† for detailed equations) as opposed to the <sup>15</sup>N isotropic chemical shifts alone. The additional <sup>14</sup>N second-order isotropic quadrupolar shifts will result in the spreading of the nitrogen signals over a larger ppm range (hundreds of ppm) and can thus disperse very similar or even overlapping <sup>15</sup>N NMR signals. This could prove to be crucial for the analysis of amorphous drug formulations, either when the nitrogen in both, the drug and polymer, is a part of the same functional moiety (*e.g.* amides in the previously mentioned POx/PTX formulations) or when the distribution of similar local environments as observed in amorphous forms makes the differentiation between peaks challenging. This is underlined by the fact that around 84% of unique small-molecule drugs contain a nitrogen atom.<sup>40</sup> As a result, there is a large, yet largely unexplored potential of this experiment for the analysis and increased understanding of such materials. In solid-state NMR, there are two different ways to create <sup>14</sup>N-<sup>1</sup>H HMQCs, either heteronuclear through-bond <sup>14</sup>N-<sup>1</sup>H  $J$ -couplings and residual second-order quadrupole-dipolar couplings or by <sup>14</sup>N-<sup>1</sup>H through-space dipolar couplings.<sup>33,34,41</sup> By varying the recoupling time in the latter case, it is then possible to obtain different through-space proximities with longer recoupling times probing longer range intra- and intermolecular



All NMR data including the compiled pulse sequences for the results presented in this paper can be found in the Warwick Research Archive Portal (WRAP). Magic angle and pulse calibrations were performed with the dipeptide  $\beta$ -AspAla.  $^{14}\text{N}$  chemical shifts were referenced to saturated  $\text{NH}_4\text{Cl}$  at  $-352.9$  ppm, corresponding to a primary reference of  $\text{CH}_3\text{NO}_2$  at  $0$  ppm.  $^{15}\text{N}$  experiments recorded at  $9.4$  T were referenced to unlabelled  $\alpha$ -glycine at  $-342.0$  ppm and those at  $20$  T were referenced to labelled L-histidine, which has an  $\text{NH}_3^+$  resonance at  $-333.1$  ppm. Both ways of  $^{15}\text{N}$  referencing correspond to a  $\text{CH}_3\text{NO}_2$  reference at  $0$  ppm. The  $^{13}\text{C}$  NMR data were measured with a Bruker Avance III HD  $600$  MHz spectrometer and a  $3.2$  mm double-channel probe. The duration of the  $^1\text{H}$  pulses was  $2.5$   $\mu\text{s}$  in the  $^{13}\text{C}$  CP MAS under ramped cross-polarization conditions optimized using a ramp from  $90$  to  $100$  with  $100$  increments and the  $\alpha$ -glycine sample. SPINAL-64<sup>46</sup> heteronuclear decoupling was applied during an acquisition period of  $23$  ms at a  $^1\text{H}$  rf nutation frequency of  $100$  kHz, with an optimized pulse length of  $4.9$ – $5.1$   $\mu\text{s}$ . The  $^{15}\text{N}$  data of POL-9-PTX were recorded with a Bruker Avance III spectrometer at a  $^1\text{H}$  Larmor frequency of  $850$  MHz and a duration of the  $^1\text{H}$  pulses of  $2.5$   $\mu\text{s}$  under the optimized cross-polarization conditions using histidine. The  $^{15}\text{N}$  data of crystalline PTX and POL were recorded using a Bruker Avance Neo  $400$  MHz spectrometer. In both cases,  $4$  mm rotors with a sample volume of  $80$   $\mu\text{L}$  were used. The duration of the  $^1\text{H}$  pulses was  $2.5$   $\mu\text{s}$  under the optimized cross-polarization conditions using  $\alpha$ -glycine. KBr was used for magic angle calibration.  $^{13}\text{C}$  chemical shifts were referenced to the methylene carbon of adamantane at  $38.48$  ppm.<sup>47</sup>


## Materials

## NMR

## Results

As the starting point for the spectral assignment and further investigation of the three differently loaded formulations of the poorly water-soluble cancer-drug paclitaxel containing 17 wt% (POL-2-PTX), 29 wt% (POL-4-PTX) and 47 wt% (POL-9-PTX) PTX, NMR experiments in solution were performed (Sections SI 4–SI 6, ESI<sup>†</sup>). For the assignment of PTX, data from NMR measurements in CDCl<sub>3</sub> were used. However, in aqueous polymer micellar formulations, only extremely broad, unresolved PTX signals can be observed in the <sup>1</sup>H NMR spectra (Section SI 6 with Fig. S4, ESI<sup>†</sup>), underlining the necessity to focus on solid-state NMR for the analysis of such formulations. Consequently, a set of different solid-state NMR spectra at moderate to fast MAS spinning frequencies were recorded for the formulations and as-received PTX, amorphous PTX and the neat polymer.

Both the  $^{13}\text{C}$  CP MAS NMR spectra recorded at 24 kHz (Fig. 1a) as well as powder X-ray diffraction data (Fig. S5, ESI $^\dagger$ ) confirm that the formulations are amorphous materials in agreement with a single glass transition temperature published in a previous work by one of the authors.<sup>9</sup> While crystalline paclitaxel (black spectrum) shows two sets of carbon resonances due to its two independent molecules in the asymmetric unit, amorphous PTX and an exemplary formulation (POL-9-PTX, highest PTX content) (grey and blue spectrum) show broad



**Fig. 1** (a)  $^{13}\text{C}$  CP MAS NMR spectra of the pure polymer (purple), an exemplary PTX formulation (blue), amorphous (grey) and crystalline PTX (black). Partial peak assignment is indicated. All spectra were recorded at 14.1 T and 24 kHz MAS with a contact time of 2.5 ms and the following parameters (co-added transients/recycle delay): 1756/2.5 s (purple), 2048/2.0 s (blue), 23 532/2.5 s (grey) and 1756/2.5 s (black). (b) Corresponding  $^{15}\text{N}$  CP MAS NMR spectra recorded at 9.4 T (PTX, POL) and 20 T (POL-9-PTX) using a contact time of 1 ms for POL and 2 ms for PTX. 82 796 (PTX), 16 991 (POL) and 25 600 (POL-9-PTX) co-added transients were measured with a recycle delay of 5.0 s (PTX, POL) and 3.0 s (POL-9-PTX) resulting in an overall experimental time of 4 d 19 h, 24 h and 21 h, respectively. Please note that the  $^{15}\text{N}$  NMR spectrum of crystalline PTX is the only dataset, which was recorded on a new batch of PTX containing two differently hydrated PTX phases and thus showing an additional set of signals. As PTX exists in its amorphous state in the formulations and is dissolved for preparation, the initial degree of incorporated water in the respective phase is not relevant for the analysis of the formulations.

signals as expected for materials lacking long-range order and containing a distribution of environments and thus chemical shifts. Consequently, several moieties of PTX appear as joint signals complicating a straightforward extraction of spectral changes upon incorporation into the micelles in a similar way as recently shown for related polymer micelles loaded with curcumin as the model compound.<sup>19</sup> Moreover, additional complications arise from the overlap between PTX and POL resonances as can be seen from a comparison with the  $^{13}\text{C}$  CP MAS spectrum of the neat polymer (purple spectrum). Note, that for better comparability, the neat polymer was heated above its glass transition temperature  $T_g$  ( $T_g(\text{POL}) = 56^\circ\text{C}$ )<sup>48</sup> prior to the measurement to account for similar conditions during the preparation of the formulations. Consequently, the analysis of the carbon spectra and extraction of reliable structural information is difficult and requires complementary tools such as quantum chemical calculations, which will be explored in detail in a separate work.

Ideally, for a thorough characterization of the formulations we would like to use the direct information from  $^1\text{H}$  NMR data at fast MAS due to the high sensitivity of the  $^1\text{H}$  nucleus to intermolecular interactions, through space proximities and subtle changes in the packing arrangement. The  $^1\text{H}$  MAS NMR spectra recorded at 60 kHz are shown as external projection in the corresponding  $^{14}\text{N}$ - $^1\text{H}$  HMQC spectra in Fig. 2 and 3. As can be seen, the spectra are dominated by severe signal overlaps despite being recorded in a high magnetic field (20 T). Consequently, spectral simplification is required and the presence of only one secondary amide moiety per PTX enables this through N-H correlation experiments. Additionally, the tertiary amide functional groups in the polymer


could also serve as another source of structural information. Therefore,  $^{15}\text{N}$  CP MAS spectra were recorded and it is very interesting to note that the secondary amides from PTX ( $Z' = 2$ , crystal structure in SI 10, ESI†) and the tertiary amide of the polymer show resonance in the same spectral region (Fig. 1b). For the POL-9-PTX formulation, only a very broad resonance can be observed, again in a similar spectral region and thus showing the need for a dispersion of the signals as expected for the quadrupolar  $^{14}\text{N}$  nucleus.

$^{14}\text{N}$ - $^1\text{H}$  HMQC MAS NMR experiments are particularly powerful in this context, require only a small amount of sample and have been applied to study a variety of crystalline compounds.<sup>33,35</sup> Therefore, in the first step,  $^{14}\text{N}$ - $^1\text{H}$  HMQC spectra of the pure compound were recorded to identify their nitrogen shifts and ensure sufficient separation of the PTX and polymer nitrogen environments.

### Pure compounds

Fig. 2 shows  $^{14}\text{N}$ - $^1\text{H}$  HMQC spectra of crystalline PTX as well as of the pure polymer. For PTX, correlations with two different dipolar recoupling times were measured and are illustrated in black (133.3  $\mu\text{s}$ ) and green (400.0  $\mu\text{s}$ ). Similar recoupling times in the  $^{14}\text{N}$ - $^1\text{H}$  HMQC experiments were first tested for the dipeptide  $\beta$ -AspAla showing the expected appearance of longer range through space proximities upon increase of the recoupling time (see Fig. S6 (ESI†) and Tatton *et al.*<sup>28</sup>). For PTX, using short recoupling times and in agreement with the X-ray diffraction data, two cross peaks at 8.8/−30 ppm ( $\text{N}^{\text{A}}$ ) and 8.0/−18 ppm ( $\text{N}^{\text{B}}$ ) were observed for the NH groups of the two individual molecules of PTX in the asymmetric unit. Consequently, one can deduce






**Fig. 2** (a)  $^{14}\text{N}$ - $^1\text{H}$  HMQC spectra (20 T, 60 kHz) of crystalline PTX recorded with recoupling times of  $\tau_{\text{RCPL}} = 133.3 \mu\text{s}$ , acquired with 34  $t_1$  FIDs (black) and  $\tau_{\text{RCPL}} = 400.0 \mu\text{s}$ , acquired with 98  $t_1$  FIDs (green, as an inset), each with 128 co-added transients for a recycle delay of 1.5 s (experimental time: 1 h 50 min (black) and 5 h 18 min (green)). The base contour level is at 30% (black), 34% (green) and 22% of the maximum peak height. The red bars point out the corresponding peaks in the  $^1\text{D}$   $^1\text{H}$  one-pulse MAS NMR spectrum. (b) The  $^{14}\text{N}$ - $^1\text{H}$  HMQC spectrum of POL recorded with a recoupling time  $\tau_{\text{RCPL}} = 266.7 \mu\text{s}$ , 16  $t_1$  FIDs, each with 2048 co-added transients for a recycle delay of 1.3 s (experimental time 1 h 30 min), and corresponding assignment of the N- $\cdots$ H cross-peaks. The skyline projections correspond to  $^1\text{H}$  one-pulse MAS NMR spectra of each sample.

that  $\text{N}^{\text{A}}$  is involved in stronger hydrogen bonding than  $\text{N}^{\text{B}}$ . This agrees with a longer direct  $\text{N}^{\text{A}} \cdots \text{H}$  distance, *e.g.* stronger involvement in hydrogen bonding, as observed in the crystal structure reported by Vella-Zarb *et al* for PTX anhydrate (Table 1, CSD code: RIGLEA).<sup>25</sup>

With a three times higher recoupling time, two changes of the previous peaks for shorter recoupling time are noticeable: the intensity of the peak corresponding to nitrogen  $\text{N}^{\text{B}}$  increases compared to  $\text{N}^{\text{A}}$ , the peak seems elongated in the  $^{14}\text{N}$  dimension and two additional signals at around 5 ppm of proton chemical shift are observed. Small changes in the  $^{14}\text{N}$  shift ( $\text{N}^{\text{A}}$ : -30 to -40 ppm and  $\text{N}^{\text{B}}$ : -18 to -22 ppm) could be caused by a temperature increase in the sample caused by the longer irradiation through the recoupling blocks. As the recoupling time for the dipolar couplings increases, it is possible to sample longer N- $\cdots$ H distances through space. Consequently, the relative signal increase of  $\text{N}^{\text{B}}$  can be explained by proximities of  $\text{N}^{\text{A}}$  and  $\text{N}^{\text{B}}$  to aromatic protons, whose chemical shifts are overlapping with those of the proton attached to  $\text{N}^{\text{B}}$  in agreement with an increased  $^1\text{H}$  intensity at 8.0 ppm. To confirm this, intra- and intermolecular close  $^{14}\text{N} \cdots ^1\text{H}$  distances of up to 3 Å were extracted from PTX crystallographic data (Table 1, CSD code: RIGLEA). The threshold was chosen based on the maximum distances previously observed for comparable recoupling times.<sup>28</sup> Indeed, the shortest contact for N1 is found for the *ortho*-CH of the adjacent aromatic ring denoted as carbon 32 in the ESI†. In turn, the shortest NH distances for the second molecule were observed for the adjacent aliphatic CH and OH units (C3'H and C2'OH), which could explain the cross peaks observed at 5.0 and 5.3 ppm in the  $^1\text{H}$  dimension. However, the distances in Table 1 can only serve as indication and should be viewed carefully as these distances originate from PXRD data at high temperature (360 K). Furthermore, it is still not fully understood why, for some small molecule compounds, fewer contacts than expected are observed.<sup>42</sup> In the  $^{14}\text{N}$ - $^1\text{H}$  HMQC MAS NMR spectrum of the neat polymer (Fig. 2b) recorded at an intermediate recoupling time of 266.7  $\mu\text{s}$  (no direct NH), three cross peaks at a  $^{14}\text{N}$  shift of 67 ppm and distinct  $^1\text{H}$  chemical shifts are observed. Peak A results from the through space proximity of the nitrogen with the methyl groups C14 and C21 as well as the methylene group C16. Peak B indicates proximity to the CH<sub>2</sub> unit of the hydrophobic butyl sidechain (C17/18) and peak C is related to N- $\cdots$ H proximities including the methyl protons of the butyl sidechain (C19). Please note, that due to micelle formation, these cross-peaks are most likely of intermolecular nature. Interestingly, no cross peaks between N and the polymer backbone CH<sub>2</sub> groups could be observed. For a more detailed discussion on the signal assignment and appearance of specific cross-peaks, the reader is referred to the SI 4 (ESI†). Overall, the comparison of the two individual components shows that their nitrogen environments can be clearly distinguished in the  $^{14}\text{N}$ - $^1\text{H}$  HMQC experiment, which is not the case if  $^{15}\text{N}$  chemical shifts are observed (see Fig. 2b and Table 1). Additionally, distinct correlation peaks are obtained in the 2D NMR spectra significantly reducing the spectral complexity in the  $^1\text{H}$  dimension. This is an essential prerequisite for the following analysis of the formulations. While revealing valuable intra- and





**Fig. 3**  $^{14}\text{N}$ - $^1\text{H}$  HMQC spectra (20 T, 60 kHz) of (a) the POL-2-PTX formulation recorded with recoupling times  $\tau_{\text{RCPL}} = 133.3 \mu\text{s}$  (light blue) and  $\tau_{\text{RCPL}} = 266.7 \mu\text{s}$  (green), both acquired with 18  $t_1$  FIDs, each with 1024 co-added transients for a recycle delay of 1.3 s (experimental time: 6 h 46 min), including tentative assignment of  $\text{N}\cdots\text{H}$  proximities. The base contour level is at 34% of the maximum peak height. The internal overlaid  $^1\text{H}$  spectrum (purple) is that for POL. (b) The POL-4-PTX formulation recorded with a recoupling time  $\tau_{\text{RCPL}} = 266.7 \mu\text{s}$  (dark blue), and 12  $t_1$  FIDs each acquired with 1200 co-added transients for a recycle delay of 1.3 s (experimental time: 5 h 17 min). The internal overlaid spectrum (orange) is the extracted slice at a  $^{14}\text{N}$  shift of  $-132 \text{ ppm}$ . The base contour level is at 22% of the maximum peak height. The skyline projections represent the  $^1\text{H}$  one-pulse MAS NMR spectra of each sample.

**Table 1** Closest  $\text{N}\cdots\text{H}$  proximities<sup>a</sup> as extracted from PTX crystallographic data (CSD code: RIGLEA, PXRD data at 363.0 K)<sup>24</sup>

$\text{N}\cdots\text{H}$ proximity	Distance/Å
$\text{N}^{\text{A}}\text{H}$	0.979
$\text{N}^{\text{A}}\text{H}-\text{C}^{\text{A}32}\text{H}$	1.642
$\text{N}^{\text{A}}\text{H}-\text{C}^{\text{A}3'}\text{H}$	2.031
$\text{N}^{\text{A}}\text{H}-(\text{C}^{\text{A}2'}\text{H}-\text{OH})$	2.601
$\text{N}^{\text{A}}\text{H}-\text{C}^{\text{A}38}\text{H}$	2.887
$\text{N}^{\text{B}}\text{H}$	0.846
$\text{N}^{\text{B}}\text{H}-\text{C}^{\text{B}3'}\text{H}$	1.915
$\text{N}^{\text{B}}\text{H}-(\text{C}^{\text{B}2'}\text{H}-\text{OH})$	2.486
$\text{N}^{\text{B}}\text{H}-\text{C}^{\text{A}31}\text{H}_3$	2.787
$\text{N}^{\text{B}}\text{H}-\text{C}^{\text{B}32}\text{H}$	2.965

<sup>a</sup> Intermolecular proximities in italics. Mobile aromatic units are shown in bold. N and H atoms of interest are underlined.

intermolecular  $\text{N}\cdots\text{H}$  proximities, the HMQC experiments also feature considerably shorter experimental times (1 h 30 min up to 7 h) than the 1D  $^{15}\text{N}$  experiments (min 1 d), while also using significantly lower samples volumes (80  $\mu\text{L}$  in 4 mm rotors *vs.* 1.5  $\mu\text{L}$  in 1.3 mm rotors).

### Formulations

Fig. 3 depicts  $^{14}\text{N}$ - $^1\text{H}$  HMQC MAS NMR spectra of two differently loaded formulations (a) POL-2-PTX and (b) POL-4-PTX, both with their  $^1\text{H}$  MAS NMR spectra as external projections. In both spectra, two different  $^{14}\text{N}$  shifts, one at positive and one at negative ppm values, can be observed with the cross-peaks at negative values being more intense than the corresponding signals at positive  $^{14}\text{N}$  shifts. For the formulation with the lowest PTX loading (POL-2-PTX), all cross-peaks appear in the aliphatic  $^1\text{H}$  chemical shift region below 2 ppm, the extracted  $^{14}\text{N}$  shifts are  $-50$  and  $+128 \text{ ppm}$ . With increasing PTX loading of the formulations, both nitrogen shifts decrease to  $-132$  and  $+103 \text{ ppm}$ . Taking a closer look at the characteristic peaks for the medium loaded formulation, POL-4-PTX, cross-peaks are again observed in the aliphatic  $^1\text{H}$  region and additional cross peaks at a proton chemical shift of 4.5 and 8.0 ppm can be clearly distinguished alongside weaker signals at 5.6 and 6.4 ppm (orange slice). The high  $^1\text{H}$  chemical shift values are indicative of aromatic or hydrogen bonding environments. Interestingly, for this sample and in contrast to POL-2-PTX, multiple, defined  $^{14}\text{N}$  shifts between  $-135$  and  $-110 \text{ ppm}$  are observed.

To complete the set of differently loaded formulations and subsequently extract trends, which might reveal information on changes of the local environment in the formulations,  $^{14}\text{N}$ - $^1\text{H}$  HMQC data were also recorded for POL-9-PTX, a formulation with almost 50 wt% PTX loading (Fig. 4). Due to hardware problems (untraceable spikes appearing in the FIDs) and despite multiple measurement attempts at different times (August 2019 and November 2019), unfortunately no artefact free dataset could be recorded. Therefore, the  $^{14}\text{N}$ - $^1\text{H}$  HMQC spectrum for POL-9-PTX was generated by addition of the individual FIDs of several datasets with different recoupling times ranging from 133.3 to 533.3  $\mu\text{s}$ . The averaging of different datasets comes at the cost of resolution, particularly in the  $^1\text{H}$  dimension. Consequently, only the  $^{14}\text{N}$  shifts and relative signal intensities






Fig. 4  $^{14}\text{N}$ - $^1\text{H}$  HMQC (20 T, 60 kHz) spectrum of the POL-9-PTX formulation. The base contour level is at 34% of the maximum peak height. The skyline projection shows the corresponding  $^1\text{H}$  one-pulse MAS NMR spectrum. Due to problems with the probe, the recorded datasets for this sample contained artefacts. Therefore, multiple datasets with different recoupling times varying from  $\tau_{\text{RCPL}} = 133.3 \mu\text{s}$  to  $\tau_{\text{RCPL}} = 533.3 \mu\text{s}$  (recycle delay of 1.0 s) were added to obtain the spectrum. Consequently, only the  $^{14}\text{N}$  shift values were used for comparison.

are discussed in the course of this work. The values  $-273$  and  $+86$  ppm are consistent with the observed decrease in shift upon increasing the PTX loading.

All extracted values for both the individual components and the formulations are summarized in Table 1, which also includes selected  $^{15}\text{N}$  chemical shifts extracted from the spectra shown in Fig. 1b. This first examination of the  $^{14}\text{N}$ - $^1\text{H}$  HMQC spectra of the three formulations is very promising as two different nitrogen environments can be distinguished for each sample, which follow a clear trend of decreasing  $^{14}\text{N}$  shift with increasing loading. An analogous trend could not have been observed based on  $^{15}\text{N}$  solid-state NMR data as indicated by the similarity of the extracted chemical shifts for the polymer and the formulation with highest PTX loading.

## Discussion

In the next step, the observed trends and key  $\text{N} \cdots \text{H}$  proximities have to be analyzed and transformed into chemical knowledge to improve our understanding of the studied formulations on the molecular level. To do so,  $^{14}\text{N}$  quadrupolar shifts  $\delta Q_{\text{iso}}$  (for detailed explanation see SI 1, ESI†) can be estimated based on the hypothesis that the  $^{14}\text{N}$  isotropic chemical shift are approximately identical to the  $^{15}\text{N}$  isotropic chemical shifts.

Table 2 Overview of the experimentally determined  $^{14}\text{N}$  shift values, selected  $^{15}\text{N}$  chemical shifts and the resulting, calculated  $^{14}\text{N}$  quadrupolar shifts

	$\delta^{14}\text{N}$ neg. ppm	$\delta^{14}\text{N}$ pos. ppm	$\delta^{15}\text{N}/$ ppm	$\delta Q_{\text{iso}}/$ ppm
Crystalline PTX	$-18/-30$	—	$-246/-260$	228/230
Amorphous PTX	$-55/-70$	—	—	200 <sup>a</sup> /185 <sup>a</sup>
Pure POL	—	+67	$-258/-263$	325/330
POL-2-PTX	$-50$	+128	—	205/383
POL-4-PTX	$-132$	+103	—	123/358
POL-9-PTX	$-273$	+86	$-263$	$-10/349$

<sup>a</sup> Calculated as  $\delta(^{14}\text{N}) - \delta(^{15}\text{N})$  using an average  $^{15}\text{N}$  chemical shift of  $-255$  ppm.

Consequently, the quadrupolar shift values can be determined by subtraction giving the values indicated in Table 2. If available, the corresponding experimental  $^{15}\text{N}$  chemical shifts were used. For all formulations, a  $^{15}\text{N}$  value of  $-263$  ppm was used. For as-received and amorphous PTX, relatively large  $^{14}\text{N}$  quadrupolar shifts of around 200 ppm are observed, which is indicative of the more asymmetric environment of the  $^{14}\text{N}$  within a secondary amide. The data for the pure polymer, which only contains tertiary amides, reveals even larger quadrupolar shifts of above 300 ppm. Interestingly, the resulting quadrupolar shifts of all three formulations remain similarly sized for the  $^{14}\text{N}$  signal observed at positive ppm values, which is further comparable in size to the neat polymer, while the quadrupolar shift of the  $^{14}\text{N}$  signal at negative ppm values is decreasing with increasing PTX loading. This indicates a more symmetric N-environment at higher PTX loadings with the very low experimentally determined quadrupolar shift for POL-9-PTX verging on that of an almost tetrahedral arrangement. How can this major change in symmetry be explained?

To approach this question, it is first necessary to discuss the assignment of the two  $^{14}\text{N}$  shifts observed for the formulations. At a first glance, the 2D NMR data of individual components, PTX and POL, showing positive as well as negative  $^{14}\text{N}$  shifts in a similar region to the two shifts found for the formulations, suggest an analogous assignment for these two environments. However, for the lowest PTX loading (POL-2-PTX), only 2.2% of all nitrogen atoms in the sample belong to PTX, while the majority of nitrogen atoms originate from the polymer. Furthermore, no intramolecular NH contact like the one found in pure PTX is observed in the formulation with the lowest loading. For POL-4-PTX and POL-9-PTX, 4.5% and 10.1% of the nitrogen atoms in the sample are found in PTX molecules. With this in mind, we hypothesize that both  $^{14}\text{N}$  signals arise from polymeric environments, the more intense signal at negative ppm values being assigned to amides in close proximity to PTX molecules acting as hydrogen bond donors and the signal at positive ppm resulting from the remaining amide fragments of the triblock copolymer. This agrees with the signal at negative ppm gaining intensity with respect to the other  $^{14}\text{N}$  environment upon increasing the PTX loading. The tertiary amides can act as hydrogen bond acceptors and in doing so the nitrogen environment gets more symmetric, which results in a smaller quadrupolar shift. For example, this has been observed for





small di- and oligopeptides, where nitrogen atoms with three bonds and thus a less tetrahedral environment due to the lone pair showed high second order quadrupolar shifts, while in the case of zwitterionic  $\text{RNH}_3^+$  environments (proton transfer being the extreme case of hydrogen bonding interaction) only a very small shift was observed.<sup>28,36</sup> Finally, the formulation with medium loading gives further insights into the drug-polymer assembly. While the amount of PTX in POL-2-PTX is very low and the  $^{14}\text{N}$ - $^1\text{H}$  HMQC is dominated by polymer-polymer contacts, distinct NH contacts with protons resonating above 4 ppm can be observed for the  $^{14}\text{N}$  at negative ppm values in the  $^{14}\text{N}$ - $^1\text{H}$  HMQC spectrum of POL-4-PTX. All these cross-peaks can only arise from intermolecular drug-polymer contacts as the polymer itself does not show any proton chemical shifts above 3.9 ppm. Specifically, cross-peaks at 8.0, 6.4, 5.6 and 4.5 ppm were observed. In comparison with the NMR spectroscopy data of pure PTX, the cross-peak above 8 ppm can be assigned to a tertiary amide  $\cdots\text{HN}$  proximity, while the remaining cross-peaks could result from interactions between the amide and different OH environments within PTX. Additionally, the spectrum of POL-4-PTX contains diverse  $^{14}\text{N}$  shifts. As this feature is neither observed for the neat polymer nor for POL-2-PTX, where only one distinct  $^{14}\text{N}$  shift is observed, it must be attributed to the increasing amounts of PTX present in the nanoparticles. Taking a closer look at the known structures of PTX with increasing number of water molecules could help to shed light on this: in the three structures published by Vella-Zarb *et al.*,<sup>25</sup> the NH group and all three OH-groups act as hydrogen bond donors to adjacent PTX and water molecules suggesting that this would also be the case for PTX-polymer interactions. The polymer amide nitrogen atoms serve as the hydrogen bond acceptor. From the previous work on similar polymers with PTX and curcumin,<sup>11,19,49</sup> we know that, with increasing loading, the poorly water-soluble guest not only interacts with the more hydrophobic pBuOx polymer block, but also with the hydrophilic pMeOx units. The critical drug loading, where the corona become significantly involved appears to be between 20–30 wt%. This would explain, why we don't observe this for the lowest loading. Furthermore, steric hindrance between the different PTX interaction sites and pBuOx *vs.* pMeOx would result in a differing deviation from ideal hydrogen bond geometries with respect to bond distances and bond angles,

while, overall, the hydrogen bonded amide nitrogen atom would be found in a more symmetric, closer to tetrahedral environment as compared to the planar,  $\text{sp}^2$  hybridized initial environment (Fig. 5).

This underlines that  $^{14}\text{N}$ - $^1\text{H}$  HMQC NMR experiments are not only spreading out the nitrogen NMR signals but are also a promising tool to study the local symmetry and detect structure determining intermolecular interactions in complex and amorphous drug-polymer formulations. To further improve the structural understanding, these finding now need to be complemented by additional experimental and detailed theoretical data.

## Summary and conclusion

In this work, we investigated amorphous polymer formulations of the anti-cancer drug paclitaxel. For both the polymer and the drug,  $^{15}\text{N}$  chemical shifts in a very similar spectral area were observed, which results in significant signal overlap, especially if amorphous samples are involved. Consequently, for this set of samples, it is very challenging to distinguish small changes in the spectra, extract trends and information on the interactions of specific moieties, *e.g.* upon increasing the drug loading of the formulations. Therefore,  $^{14}\text{N}$ - $^1\text{H}$  HMQC NMR experiments were investigated as a valuable tool to disperse signals due to the additional quadrupolar shift when observing  $^{14}\text{N}$  instead of  $^{15}\text{N}$ . We could show that such experiments can be a rich source of information for the characterization of molecular interactions in amorphous polymer-drug nanoformulations, which is otherwise difficult to obtain. On the one hand, the  $^{14}\text{N}$  signals were dispersed so that two distinct  $^{14}\text{N}$  environments could be extracted for each formulation and the complexity in the  $^1\text{H}$  dimension could be reduced through this correlation experiment. Thus, it was possible to extract trends for different drug loadings, distinguish cross peaks and extract quadrupolar parameters, which altogether enabled us to learn about the local amide environment, its symmetry and potential interacting motifs between the drug and polymer in such complex amorphous mixtures. Overall, the  $^{14}\text{N}$ - $^1\text{H}$  HMQC NMR experiments have great potential for the analysis of disordered and amorphous drug-delivery systems, which can be transferred and should be explored for other systems. Here, the  $^{14}\text{N}$ - $^1\text{H}$  experiments could also be insightful for nanocrystalline materials such as those found in Pluronic F127 stabilized PTX nanocrystals.<sup>50</sup>

It is of importance to stress for the broader applicability of the  $^{14}\text{N}$ - $^1\text{H}$  HMQC experiment that while fast MAS ( $\geq 60$  kHz) is essential for these experiments due to the narrowing of the  $^1\text{H}$  signal and most importantly increased  $^1\text{H}$  coherence lifetimes,<sup>28</sup> the experiments are not limited to high magnetic fields. With the second order quadrupolar interaction being inversely proportional to the strength of the magnetic field, narrower lines are obtained at higher fields. However, a larger dispersion of the signals is counterbalancing this at lower magnetic field strengths. In fact, many examples for the application of this and other  $^{14}\text{N}$  based experiments in the literature show their feasibility at 400–600 MHz.<sup>30,34,36,51</sup>




Fig. 5 Summary of the hydrogen bond donating moieties in PTX alongside nitrogen atom geometries and environments for initial tertiary amides and hydrogen bonded functionalities.





- 37 A. S. Tatton, T. N. Pham, F. G. Vogt, D. Iuga, A. J. Edwards and S. P. Brown, *CrystEngComm*, 2012, **14**, 2654–2659.
- 38 A. S. Tatton, H. Blade, S. P. Brown, P. Hodgkinson, L. P. Hughes and S. O. N. Lill, *et al.*, *Cryst. Growth Des.*, 2018, **18**, 3339–3351.
- 39 D. Bernasconi, S. Bordignon, F. Rossi, E. Priola, C. Nervi and R. Gobetto, *et al.*, *Cryst. Growth Des.*, 2020, **20**, 906–915.
- 40 E. Vitaku, D. T. Smith and J. T. Njardarson, *J. Med. Chem.*, 2014, **57**, 10257–10274.
- 41 S. Cavadini, A. Abraham and G. Bodenhausen, *Chem. Phys. Lett.*, 2007, **445**, 1–5.
- 42 S. P. Brown, *High-resolution 1H 2D Magic-angle Spinning Techniques for Organic Solids. Modern Methods in Solid-state NMR, A Practitioner's Guide: The Royal Society of Chemistry*, 2018, ch. 2, pp. 39–74.
- 43 M. M. Lübtow, L. Hahn, M. S. Haider and R. Luxenhofer, *J. Am. Chem. Soc.*, 2017, **139**, 10980–10983.
- 44 T. Oas, R. Griffin and M. Levitt, *J. Chem. Phys.*, 1988, **89**, 692–695.
- 45 P. Costa, J. Gross, M. Hong and R. Griffin, *Chem. Phys. Lett.*, 1997, **280**, 95–103.
- 46 B. Fung, A. Khitrin and K. Ermolaev, *J. Magn. Reson.*, 2000, **142**, 97–101.
- 47 C. R. Morcombe and K. W. Zilm, *J. Magn. Reson.*, 2003, **162**, 479–486.
- 48 M. M. Lübtow, M. S. Haider, M. Kirsch, S. Klisch and R. Luxenhofer, *Biomacromolecules*, 2019, **20**, 3041–3056.
- 49 M. S. Haider, M. M. Lübtow, S. Endres, V. Aseyev, A.-C. Pöppler and R. Luxenhofer, *ACS Appl. Mater. Interfaces*, 2020, **12**, 24531–24543.
- 50 J. Deng, L. Huang and F. Liu, *Int. J. Pharm.*, 2010, **390**, 242–249.
- 51 D. Carnevale, X. Ji and G. Bodenhausen, *J. Chem. Phys.*, 2017, **147**, 184201.

